




Medela Clario Toni Home Care pump



1





2






3

min. -70 mmHg / -9 kPa
(+/-20 mmHg / +/-2.6 kPa)
med. -115 mmHg / -15 kPa
(+/-30 mmHg / +/-4 kPa)
max. -600 mmHg / -80 kPa
(-90 mmHg / -12 kPa)





* Measured at 0m, atmospheric pressure 1013.25 hPa. Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).


4



5



6



This Quick Card does not replace the instructions for use REF 200.8400/200.8401



Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medelahealthcare.com

CE 0123

Precious life – Progressive care